

# Thermocycling Run on TEC A


This is a report automatically generated after a run on TEC A. The contents of this report include when the run started and ended, which protocol was used, plots for the expected temperature, actual temperature (where actual here means the temperature measured by the internal temperature sensor for TEC A), the sink temperature, and when the fans were on for TEC A. Extra information regarding the cartridge used, the pressure applied to said cartridge, the elastomer/bergquist used if at all, images before/during/after, and positions for the clamp and tray are not included within this report.

## Protocol Status


This section includes the protocol statuses for all TECs at the end of this run.

	TEC A	TEC B	TEC C	TEC D
State	Connected	Not Connected	Not Connected	Not Connected
Protocol	Sample MW dPCR Protocol	?	?	?
Protocol Running	No	?	?	?
Estimated Time Left	00:00:05	?	?	?
IO	Idle	?	?	?
Step	5	?	?	?
Cycle Number	3/3	?	?	?
Temp (°C)	104.0049	?	?	?
Target Temp (°C)	30.0000	?	?	?
Sink Temp (°C)	35.8616	?	?	?
Fan RPM	10275.2021	?	?	?

## TEC A Object and Target Temperatures



## TEC A Sink Temperatures



**TEC A Fan Speeds**

